

DUAL WAVELENGTH FOR ULTIMATE CONVENIENCE




755nm
1064nm

What our YAG do for you? **Versatility!**

Our 755 nm and 1064 nm wavelength technology goes beyond hair removal to also let you treat wrinkles, benign pigmented lesions and vascular lesions. That flexibility can lead to new revenue and boost your ROI.




WORKING PRINCIPLE




- During Treatment

Non-invasive laser energy is absorbed by hemoglobin, and converted to heat energy.



- Immediately After Treatment

Blood coagulates and the vessel wall collapses.



- After Treatment

Over time, the vessel will be reabsorbed by the body and disappear.

TREATMENTS



Hair Removal



Hemangioma
Removal



Blood Vein
Removal



Vascular Removal



Skin Rejuvenation



Nail ringworm
Treatment



HAIR REMOVAL

1064NM ND:YAG LONG PULSE

Ferrari in hair removal

THE MOST EFFECTIVE
HAIR REMOVAL DEVICE

