



eXodus
Lasers


HIFEM+EMS+RF

Three therapies in a single treatment.

BUILD MUSCLE

SCULPT YOUR BODY

SKIN TIGHTENING



Renu Sculpt Plus

Renu Sculpt Plus builds off the legacy of its predecessor, Emsculpt, by simultaneously emitting both Radio Frequency and High-Intensity Electromagnetic Energies.

The end result is more fat reduction and muscle growth than any single gold standard procedure, for less time and less money.



The Renu Sculpt Plus System Highlights

Why is Renu the ONE in body shaping?

Renu is the specially and only non-invasive body shaping procedure that combines 2 energies for simultaneous fat elimination and muscle building in a 30-minute session.

What areas can be treated with Renu Sculpt Plus

Renu can be applied over the abdomen, buttocks, thighs, arms & calves.

What does Renu therapy feel like?

You can feel intense muscle contractions together with a heating sensation in the treated area which is comparable to hot stone massage.

How does Renu compare to liposuction and non-invasive device alternatives?

Renu does not compare to liposuction. As FE60 is a completely non-invasive treatment. It does not require surgery, needles, anesthesia, or any downtime. Organically combines HIFEM and RF to provide three therapeutic effects: muscle building, fat reduction and skin tightening during one treatment.

Workout Smarter Not Harder

Non-essential fat provides greater resistance to athletic motion thereby forcing one to increase the muscle force of contractions per given workout. This takes time, energy and a lot of patience, until one discovers Renu Sculpt Plus



High-Intensity Focused Electro-Magnetic Technology

HIFEM technology, the full name is High-Intensity Focused Electro-Magnetic, it induce short bursts of powerful muscle contractions that lead to increased muscle density, decreased volume, better definition, and improved tone. Also the energy will eliminates the fat cells, so not only for muscle building, but also for fat reduction. Clinical research shows that after a course treatment, it can increase 16% muscle and get 19% fat reduction.



Electro Muscle Stimulation

EMS is a type of electrotherapy stimulates a muscle contraction using electrical impulses, this stimulation creates muscle contractions that can be quick and frequent, fast with long pauses, so get good results on muscle building as well.



Radio Frequency

RF is a nonsurgical method of tightening your skin. The procedure involves using energy waves to heat the deep layer of your skin known as your dermis. This heat stimulates the production of collagen so that get skin firming and tightening results.



RF+EMS



HIFEM

30% Fat burn

HI-FEM technology combines micro EMS+ bipolar RF technology, simultaneously send out high-frequency contraction signals, which directly destroys fat cells and makes them die naturally. At the same time, Due to the radiofrequency heating, the fat cells are permanently damaged and cleared from the body. Maintains warmth and comfort during treatment, creating beautiful body lines.

25% Muscle gain

HIFEM combines with EMS and RF contracts the muscle fibers in this area at a strength and improve the blood circulation of the body .In 30 minutes, increase in the number and growth of muscle fibers and cells and the muscle building effect will be seen immediately. 4 sessions equals 12-16 weeks of exercise.



Combines 3 energies for simultaneous fat elimination and muscle building in a 30-minute session.

Multiple working modes, Provide diverse Therapeutic experience
Two HI-FEM + EMS&RF working modes

- ① Combined mode, three technologies work at the same time.
- ② Alternate mode: HI-FEM first, then EMS&RF, cycle work



Before



After



Big Handle × 2

You begin to feel tangible results right after the treatment. Positive results are usually reported two to four weeks after the last session and continue to improve for several weeks following the treatments.



Small Handle × 2

Renu Sculpt Plus treatment does not require any pre-treatment preparations. You can follow your regular lifestyle right after the treatment.



Before



After



Pelvic Chair Cushion

- 1- The active regeneration of collagen in pelvic floor tissues, strengthens pelvic floor muscles, and lifts buttocks
- 2- Improve privacy issues, such as urinary incontinence, penetration, privacy issues
- 3- Achieve vagina tightening.



Before



After



Product Name	FE60
Technology	HIFEM+RF+EMS
Handle Piece	2 Big Size & 2 Small Size & Pelvic Chair Cushion
Magnetic Vibration Intensity	1-13 Tesla
RF Frequency	20MHZ
HIFEM Frequency	1-135HZ
Treatment Mode	Combined Mode& Alternate Mode
Input Voltage	AC110V-220V
Output Power	5000W
Pulse Duration	350µs
Fuse	20A
Cooling System	Air Cooling
Size of Flight Shipping Case	69*57*112cm
Gross /Net Weight	90/71KG
Working Mode	Automatic Mode and Manual Mode

shell materials are made of first-class

Platform Benefits

* 3 technologies in 1 — HIFEM+EMS+RF.

The shell materials are made of first-class

ABS. The machine only occupies small area and space with compact structure.

* Dual effects — achieving muscle building and fat reducing at the same time.

* 5 handles optional — for independent and synchronized work on different body parts.


* Newly-added cushion — easily reducing fat on buttocks while sitting.

* 10-inch touch screen — simple to use with clear guidance of intellectualized design.

* Time & effort saving — lying / sitting down for 30min = 3000 muscle contractions = 3000 belly rolls / squats.

* Remarkable results — a course of treatment can effectively increase muscle by 25% and reduce fat by 30% at the same time.

* Safe & comfortable — non-invasive treatment with no radiation, pain or side effects.



Exodus
Lasers

