



Platinum Milk light

A Hair Removal Device that Can Rejuvenate the Skin



WORKING PRINCIPLE

Platinum Milk light



Platinum Milk Light is the first comprehensive skin rejuvenation machine in China that integrates milk light skin beautification + hair removal, which can repair various skin problems of customers. The superimposing action makes the skin multi-dimensionally rejuvenated, doubles the clarity, and instantly flashes the beautiful skin.

· Milk light skin beautifying function ·

The selected near-infrared precision light has a wavelength of 900nm – 1800nm and peaks at 1300nm. It can effectively reach the skin root that is 0.5 – 2mm under the skin. It can also reach the middle to deep layer of the dermis to speed up the cell metabolism., improve the skin texture, increase the skin surface refractive index, and make the skin glow from the inside out.

· Depilation function ·

The skins with different hair and color tone have different absorption rates at different laser wavelengths. The three wavelengths of 808nm, 755nm, and 1064nm are perfectly combined to precisely act on the specified target and perform high-energy and high-precision actions on the root of the hair follicles, making the hair shed automatically and truly achieving the permanent hair removal with zero pain.



7

MAJOR EFFECTS

Platinum Milk light



1 Improve skin tone

2 Shrink pores

3 Fight anti-aging

4 Hair removal

5 Lighten stains

6 Reduce redness

7 Stimulate collagen regeneration


· Milk light handle ·

(body handle, face handle: two light spots can be freely replaced)

- ⌚ Intelligent low temperature cooling, the whole process is comfortable and painless.
- ⌚ Perfect pulse, effective and full of energy.
- ⌚ Innovative sliding, uniform and full of light energy coverage.
- ⌚ Germany Heraeus special light source ensures the energy stability and safety.



NIR handle spot size:
8* 40mm / 20*60mm



Energy: 10-140J
Power: 800 / 1500 / 2000



Hair removal handle spot size:
15*15mm / 15*30mm

· Hair removal handle ·

- ⌚ China's first 365g light-weight hand tool, which provides efficient hair removal.
- ⌚ Three wavelengths cover all-round services for all skin tones.
- ⌚ The integrated sapphire crystal filter provides more accurate light output.
- ⌚ The top-grade laser component has 72 hours longer lifetime.



Exodus Lasers

Reshape high-gloss skin with the light of technology



Platinum Milk light

FEATURE HIGHLIGHTS



Whitening + hair removal dual effect

The energy is more concentrated, which improves the skin texture, skin whitening, and hair removal. The skin rejuvenation and whitening effect can be quickly achieved. The efficiency and function of the machine can also be greatly improved.



Efficient

The sliding technology, combined with the ultra-high frequency, releases the energy evenly into the treated area, ensuring more comprehensive and effective coverage of all skins.



Precise arrival

The near-infrared precise light waves can directly act on the dermis, improving the functions of the body cells and instantly making the skin healthy.



Comfortable and safe


The energy is precisely applied to the target tissue, which only promotes the absorption of excess heat generated by water molecules. Thus, other tissues of the epidermis are not affected. This makes the heating even and safe.

Platinum Milk light

TREATMENT EFFECT

EXXodus
Lasers

360° instant radiance
of the whole
body skin



Face

Neck

Limbs

Back

Hair removal

Waist

Abdomen

Buttocks

Bikini

MARKET DEMAND

1. The first machine to cover the maintenance needs of the whole body skin.
2. All-round skin beautifying experience, hair removal + flawless skin, lasting smoothness and whitening.