



EXCOUS

# Exo Prime le

A breakthrough hair removal platform



Accurate three-wavelength For all skin types



Super Power



Big spot size (15\*30mm)









# Exo Prime Ice

For all skin types)

Exo Prime Ice uses 3000W high-power precision three-wavelength treatment technology. The light energy directly acts on the hair follicle tissue of the dermis. It removes the melanin from the hair follicle tissue without damaging normal skin and sweat glands, and achieve the purpose of ultra-fast and 100% painless hair removal.

A breakthrough hair removal platform that combines ultimate functionality with unparalleled comfort. Exo Prime Ice leverages the benefits of three combined wavelengths, with an Improved patient experience and a business-oriented approach. creating a unique and result driven new solution in the world of professional hair removal.





fter

①Emit a laser beam. ②The laser beam is absorbed by hair follicle melanin. ③Different color hair has different ability to absorb light wave. ④After absorption, the laser energy will be converted into heat energy. ⑤Hair follicles and surrounding stem cells lose their activity effect and fall off. ⑥So reach hair removal effect.



3000w Super high power



Simultaneous contact cooling



Three waves diode laser (755nm+808nm+1064nm)



Support 24 hours continuous work



## Exo Prime Ice

For all skin type

#### The most popular non-invasive treatment

- Painless: The unique epidermal protection mode can effectively protect the skin from damage, it is comfortable without burning and tingling.
- Fast: Change the treatment dead zone formed by the traditional laser circular spot, which can remove hair in a large area, with good effect and fast speed.
- Safety: Directly hit the melanin in the hair follicle, the epidermis absorbs less energy, does not damage the surrounding tissues, and has no inflammation and pigmentation.
- Precise: The light energy is directly absorbed by the melanin particles in the hair and hair follicles. The light and heat effect completely destroys the hair follicles and hair shafts, so as to effectively remove hair.

#### Handle and advantage

Ultra-light hand with light texture, the latest generation of laser technology. Fast treatment of large spots, treatment time and course reduction

- High-power 3 wavelengths to meet the needs of whole body hair removal
- High frequency, fast and smooth improvement of treatment efficiency
- Upgrade the cooling system to effectively ensure 24 hours non-stop





### Why choose 3000W diode laser?

3000W diode laser from Exo Prime Ice stands for the highest technology in the hair removal industry, which can help you achieve super fast, permanent and 100% painless hair removal.









